

Facility Name: _____

Clinical Trainee Name: _____ Training Date: _____

Trainer's Name: _____

Circle One: Satisfactory Completion Further Review Needed

Training Goals	Trainer Initials
1. Describe the intended use and contraindications of the BiWaze Clear System	
2. Explain the mechanism of action for the OLE therapy delivered by BiWaze Clear	
3. Assemble the Dual Lumen Breathing Circuit, connect to AC power and switch on the control unit	
4. Demonstrate how to fill the Aerogen® Solo nebulizer and insert the Aerogen power cable	
5. Review the proper patient positioning and technique for each patient interface	
6. Describe the Auto therapy mode and Manual therapy modes as well as the benefit of therapy profiles	
7. Demonstrate how to program a Manual therapy and describe the different phases	
8. Demonstrate how to run a Manual therapy and how to adjust pressures and frequency while the therapy is being delivered to the patient	
9. Demonstrate how to program an Auto therapy for a typical treatment and how to save therapy profiles.	
10. Demonstrate how to run an Auto therapy and describe the information being displayed on the screen.	
11. Explain how Manual and Auto mode are different and the benefits of each therapy mode.	
12. Explain why you would lock the control unit and demonstrate how to lock the control unit	
13. Describe which components stay with the control until and which are disposable	
14. Describe the setup steps for placing BiWaze Clear in-line with the ventilator	
15. Describe which mode and pressure should be used with BiWaze Clear and what ventilation mode should be used when placing BiWaze Clear in-line.	
16. Describe what to do if the patient is showing signs of respiratory distress.	
17. Demonstrate what to do when the BiWaze Clear therapy is complete.	
18. Explain why BiWaze Clear should never be left in the spring-valved tee adapter if the therapy isn't running.	
19. Describe how to clean the control unit and store the Dual Lumen Breathing Circuit	

Clinician Signature: _____

Trainer Signature: _____