

BiWaze[®] Clear System – effects of oxygen bleed-in on FiO₂ delivery in vitro

Niko Kontoudios RRT , Robert DiBlasi RRT-NPS, FAARC

Seattle Children’s Hospital, Seattle, Washington, USA

INTRODUCTION

In the field of respiratory therapy, oxygen (O₂) therapy plays a critical role in the management of patients with compromised respiratory function. One essential parameter in O₂ therapy is the fraction of inspired oxygen (FiO₂), which refers to the concentration of O₂ in the inspired gas mixture. Accurate control and maintenance of FiO₂ levels are crucial in providing optimal oxygenation while minimizing the risk of O₂ toxicity or hypoxia. O₂ bleed-in is a technique that uses an auxiliary gas flow to supplement the delivered FiO₂ during oscillating lung expansion (OLE) therapy.



This bench study aimed to investigate the effects of O₂ bleed-in on FiO₂ levels during OLE therapy with BiWaze Clear (ABM Respiratory Care, USA) under various therapeutic parameters.

STUDY METHOD

We utilized the ASL 5000 spontaneously breathing lung simulator (Ingmar Medical, USA) to simulate two patient models: an adult patient weighing 70 kg and a pediatric patient weighing 20 kg. The BiWaze Clear system was configured to deliver two therapy phases: positive expiratory pressure (PEP) and high-frequency oscillation (OSC). Each phase involved the evaluation of two pressure settings, representing a range commonly used during OLE therapy.

FiO₂ measurements were obtained for each O₂ bleed-in flow setting, starting from 1 liter/minute and incrementally increasing up to 15 liters/minute. O₂ bleed-in was achieved by introducing supplemental O₂ flow into the inspiratory path of the coaxial bacterial/viral filter. A paramagnetic O₂ sensor integrated into the ASL lung model was used to continuously analyze the FiO₂ level.

RESULTS

Adult with PEP of 5 cm H ₂ O	
Bleed O ₂ (L/min)	Delivered FiO ₂ (%)
1	24.00
2	29.00
3	35.00
4	41.60
5	48.60
6	54.60
7	61.60
8	67.10
9	76.10
10	81.70
11	86.40
12	89.50
13	91.80
14	93.30
15	94.10

Adult with PEP of 15 cm H ₂ O	
Bleed O ₂ (L/min)	Delivered FiO ₂ (%)
1	24.50
2	27.50
3	31.50
4	35.00
5	38.50
6	42.00
7	45.50
8	49.00
9	52.50
10	56.00
11	59.50
12	63.00
13	66.50
14	70.00
15	73.50

Adult with OSC of 10 cm H ₂ O at 4 Hz	
Bleed O ₂ (L/min)	Delivered FiO ₂ (%)
1	24.50
2	27.50
3	31.50
4	35.00
5	40.00
6	44.00
7	48.00
8	51.50
9	55.00
10	58.50
11	62.50
12	66.00
13	69.50
14	73.00
15	76.50

Adult with OSC of 30 cm H ₂ O at 4 Hz	
Bleed O ₂ (L/min)	Delivered FiO ₂ (%)
1	23.30
2	25.50
3	28.30
4	31.40
5	34.60
6	37.80
7	40.70
8	43.80
9	46.20
10	49.10
11	51.50
12	54.70
13	57.20
14	59.30
15	62.20

Pediatric with PEP of 5 cm H ₂ O	
Bleed O ₂ (L/min)	Delivered FiO ₂ (%)
1	26.70
2	34.30
3	43.00
4	49.60
5	55.40
6	63.00
7	69.80
8	74.80
9	80.00
10	84.30
11	88.00
12	91.60
13	93.30
14	94.50
15	94.90

Pediatric with PEP of 15 cm H ₂ O	
Bleed O ₂ (L/min)	Delivered FiO ₂ (%)
1	26.00
2	32.50
3	38.00
4	44.00
5	50.00
6	56.00
7	60.00
8	66.00
9	71.00
10	76.00
11	80.00
12	84.00
13	87.50
14	90.00
15	92.50

Pediatric with OSC of 10 cm H ₂ O at 4 Hz	
Bleed O ₂ (L/min)	Delivered FiO ₂ (%)
1	26.00
2	32.50
3	40.00
4	46.00
5	53.00
6	58.00
7	64.00
8	70.00
9	75.00
10	79.00
11	84.00
12	88.00
13	91.00
14	93.50
15	95.00

Pediatric with OSC of 30 cm H ₂ O at 4 Hz	
Bleed O ₂ (L/min)	Delivered FiO ₂ (%)
1	24.50
2	26.80
3	29.90
4	33.80
5	37.80
6	41.60
7	45.10
8	48.60
9	52.30
10	57.80
11	58.10
12	60.80
13	63.90
14	66.60
15	68.70

CONCLUSION

This bench study aimed to evaluate the effects of O₂ bleed-in on FiO₂ levels with a variety of therapeutic parameters during OLE therapy with BiWaze Clear. Through a comprehensive evaluation of different therapy parameters, our analysis showed a linear correlation between O₂ bleed-in flow and delivered FiO₂ for all therapy phases and pressures. We created guidance tables with the measured values of FiO₂ for the analyzed bleed-in flow rates for both patient models during the different therapy phases and settings. The findings of this study can aid clinicians in selecting the appropriate flow to bleed-in to optimize FiO₂ delivery with BiWaze Clear.

BiWaze® and ABM Respiratory Care® are registered trademarks of ABMRC, LLC.

©2023 ABMRC, LLC. ALL RIGHTS RESERVED. PRTN-1585539264-714 Rev. 1.0 JUL-2023 ENG US

For further information about this product or a service, visit our webpage: www.abmrc.com or call customer service: 1-877-ABM-RC01 (877-226-7201)

